


# HYPOTONIC + HYDROGEN



## The Vital power “Ultra Pure Water”+“Hydrogen”

Aqua Street's H+H can offer you  
high concentrated Hydrogen Pure Water



### Hardness 0

It means Hypotonic which character has  
high absorption & penetration to human body.



### High Concentrated Hydrogen

Contains 1.40ppm Hydrogen



### Long time preservation

Packaged in QUAD Aluminum Laminated Pack + Patented  
Hydrogen Filling System (Patent No.3606466 in JAPAN)



NEW

Aqua Street

# HYPOTONIC+HYDROGEN WATER

High Concentrated Hydrogen  
0.80ppm~1.40ppm



Aqua Street

# Characteristics of Aqua Street's H+H Water



## Hardness 0

**High Rejection:** Inorganic minerals, heavy metals, chemical compounds are highly rejected by RO system.

## Hypotonic

**Quick Penetration:** Our hydrogen water has a distinctive character to penetrate into human body quickly. It has well known by providence nature, so called "Osmosis".

## High Conc.

**1.40 ppm:** Our filling system can realize to fill up hydrogen to 1.60ppm which is theoretically max. Even after heat treatment 1.40ppm is confirmed.

## Exact Sealed

**QUAD Alminum:** Completely sealed by QUAD aluminum laminated pack can offer to maintain high concentrated hydrogen 0.80ppm even 3 month after dispatched our factory.

## Preservation

**24 Hours:** Even cap opened, you can enjoy high hydrogen water for 24 hours by pushing out remaining air from pack.




## HYPOTONIC+HYDROGEN WATER

- Aqua Street has been contributing to develop RO Water Purifier for 29 years. The pure water vending machine well known is our leading-edge technology.
- We contribute to realize your healthy life in long range by collaborating with RO water purification systems and high concentrated hydrogen.

### Nutrition Information PER 100ml

● Energy .....	0	● Magnesium .....	0
● Fat .....	0	● Potassium .....	0
● Sodium .....	0	● Hardness .....	0
● Calcium .....	0	● Carbohydrate .....	0
● HYdrogen (at dispatched) .....	1.40ppm		
● ORP (at dispatched) .....	-550mV~ -650mV		


### Measuring Test Result of Hydrogen-1



#### Test Result

Confirmed high concentration until 3 month after filled up. Concentration will be reduced after 3 month.

### Measuring Test Result of Hydrogen-2



#### Test Result

Only 3% reduce of hydrogen in 16 hours later after re-capped and air rejected from pack. Consideration: Existing of high concentration of hydrogen are confirmed 24 hours later.

Sales Agency / **Kankyo Kougaku Inc.**  
2-40-11 Kami-ikedai Oota-ku Tokyo Japan  
145-0064

TEL **+81-03-6425-2400**

<http://www.aqua-street.com>